

FORCE TENSIOMETER – K6

SPECIFICATIONS



Product group specifications**K6****Force measurement**

Calibration	external weight
Calibration weight	CP11 or CP12 (optional)

Sample stage

Travel distance	85 mm
Simple platform	yes

Measurement specifications**K6****Du Noüy ring**

Results	surface tension (SFT)/interfacial tension (IFT)
Range	1 to 90 mN/m
Resolution	0.5 mN/m
Correction method	linear correction

General specifications**K6****Housing and peripherals**

Built-in bubble level	yes
-----------------------	-----

Instrument dimensions

Footprint	230 mm × 270 mm (W × D)
Height	315 mm
Weight (without accessories)	3 kg